

EES Product Catalog 2012



Ethicon
Endo-Surgery


PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES



Ethicon

Endo-Surgery

Transforming patient
care through innovation.




Ethicon Endo-Surgery 2011 Innovations




EES Generator




HARMONIC SYNERGY* Combination Hook Blade



ECHELON FLEX™ 45mm & 60mm
ENDOPATH™ Staplers



ETHICON SECURESTRAP™



ENSEAL™ G2 Super Jaw Tissue Sealer




HARMONIC FOCUS®
Long Curved Shears




ENDOPATH® XCEL™
Bladeless Trocar



LIGAMAX™ 5




LAPRA-TY® Clip Applier KA200L



ECHELON FLEX™ Powered ENDOPATH® Stapler

CONTENTS

Our Story	7
Energy	9
Access	19
Stapling	33
Ligation	55
Endo Devices	61
Gastric Band	71
Ordering	76
Index	77



Ethicon Endo-Surgery

EES is part of the Johnson & Johnson Family of Companies. We not only represent a history of innovation, but a rich tradition of patient care that begins with our parent company, Johnson & Johnson.

The mission of EES is to improve patients' lives by developing innovative technology-based solutions that deliver superior outcomes.

- We are guided by the Johnson & Johnson Credo — putting customers first
- We are a leader within J&J Surgical Care — part of the J&J Medical Devices and Diagnostics sector
- We are recognized as an innovator in the healthcare industry
- We have touched the lives of over 17 million patients worldwide over the past five years






Energy


EES ENSEAL®

ENSEAL® G2 Super Jaw Tissue Sealer




	CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
NEW	NSEALX22L	40mm Curved Jaw	12mm	22cm	6

ENSEAL® TRIO





	CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
	ETRIO314H	3mm Curved	5mm	14cm	6
	ETRIO325H	3mm Curved	5mm	25cm	6
	ETRIO335H	3mm Curved	5mm	35cm	6
	ETRIO345H	3mm Curved	5mm	45cm	6

ENSEAL® Round Tip



	CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
	NSEAL514RH	5mm Round	5mm	14cm	6
	NSEAL525RH	5mm Round	5mm	25cm	6
	NSEAL535RH	5mm Round	5mm	35cm	6
	NSEAL545RH	5mm Round	5mm	45cm	6

ENSEAL® Standard Tip



	CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
	NSEAL514H	5mm Standard	5mm	14cm	6
	NSEAL525H	5mm Standard	5mm	25cm	6
	NSEAL535H	5mm Standard	5mm	35cm	6
	NSEAL545H	5mm Standard	5mm	45cm	6

Strong on Sealing. Gentle on Tissue.

For your open procedures, the new ENSEAL® G2 Super Jaw Tissue Sealer brings a new level of confidence and consistency to sealing large vessels up to 7mm.

The next generation ENSEAL® technology delivers sealing advancements:

- **Strong:** Seals vessels with seal strength that withstands up to seven times normal systolic pressure¹
- **Complete:** Delivers high uniform compression throughout the 40mm curved tip²
- **Contained:** Provides minimal lateral thermal spread approximately 1mm, secure tissue grasping and minimal sticking



¹Tested in porcine carotid arteries. Normal systolic = 120mm Hg. Data on file.
²Compression testing conducted on inanimate model. Data on file.
© 2011 Ethicon Endo-Surgery, Inc. All Rights Reserved. DSL 10-1336.CR

For more information, contact your
Ethicon Endo-Surgery Sales Professional or call 1-800-USE-ENDO.

Many procedures are complex. This one is simple.

The new Ethicon Endo-Surgery Generator combines advanced technology, multifunctionality and touchscreen simplicity in one easy-to-operate unit that supports current and future HARMONIC® and ENSEAL® devices.

- Technology platform to support current and future HARMONIC® and ENSEAL® devices
- Touchscreen for fast and easy setup, operation and onscreen diagnostics
- High resolution display with wide viewing angles
- Software updates via USB memory stick




SHEARS

HARMONIC ACE®





CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
ACE23E	Curved Shears	5mm	23cm	HP054	6
ACE36E	Curved Shears	5mm	36cm	HP054	6
ACE45E	Curved Shears	5mm	45cm	HP054	6

HARMONIC ACE®





CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
ACE14S	Curved Shears	5mm	14cm	HP054	6

HARMONIC FOCUS®



CODE	TIP	SHAFT LENGTH	HANDPIECE	QUANTITY
FCS9	Curved Shears	9cm	HPBLUE	6
FCS17	Curved Shears	17cm	HPBLUE	6


HARMONIC WAVE®



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
WAVE18S	Straight Shears	8.5mm	18cm	HP054	6



EES HARMONIC®

HARMONIC® Coagulating Shears



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
LCSC5	Curved Shears	5mm	36cm	HP054	6


HARMONIC® Coagulating Shears



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
CS14C	Curved Shears	5mm	14cm	HP054	6
CS23C	Curved Shears	5mm	23cm	HP054	6

BLADES

HARMONIC SYNERGY®



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
SNGCB	Curved	3mm	4-9cm	HPBLUE	6
SNGHK	Hook	3mm	4-9cm	HPBLUE	6
SNGHK2	Hook	3mm	4-9cm	HPBLUE	6


EES HARMONIC®

HARMONIC® Laparoscopic Blades



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
HC325	Curved	5mm	32cm	HP054	6
HDH05	Hook	5mm	32cm	HP054	6

HARMONIC® Ball Coagulator



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
HBC05	Ball	5mm	31cm	HP054	6

HARMONIC® Blade



CODE	TIP	SHAFT DIAMETER	SHAFT LENGTH	HANDPIECE	QUANTITY
DH105	Hook	5mm	10cm	HP054	6

EES HARMONIC®

Endo-Surgery Generator



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	GEN11	EES Generator	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
	HP054	HARMONIC® Handpiece	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
	HPBLUE	HARMONIC® Handpiece (Blue)	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	HGA11	Connector: compatible with legacy HARMONIC® Handpieces	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	EGA11	Connector: compatible with legacy ENSEAL® Instruments	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	GEN11VK	Output Verification Key	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	FSW11	Footswitch compatible with EES Generator	1



	PRODUCT CODE	DESCRIPTION	QUANTITY
NEW	CRT11	Cart for EES Generator	1

ACCESSORIES

HARMONIC® Accessories

CODE	DESCRIPTION	QUANTITY
HST01	Sterilization Tray	1
HST02	Test Tip for HP054	1
TLB01	Torque Lock Blade Wrench	2
TTBLUE	Test Tip for HARMONIC® Handpiece (Blue)	1
HSA06	5mm Blade System Adapter	2
HSA08	Disposable Hand Switch Adapter (only use with HP054)	3
HAOSSL	HARMONIC ACE® 5mm Open Sleeve	6

WARRANTIES

EES Energy Generator Extended Warranty

CODE	DESCRIPTION	QUANTITY
GEN11WAR1	EES Energy Generator 1-Year Extended Warranty	1
GEN11WAR2	EES Energy Generator 2-Year Extended Warranty	1





Access

New ENDOPATH® XCEL™ with OPTIVIEW® Technology

reduces smudging during scope reinsertion more effectively than any other trocar.*

With the addition of NEW OPTIVIEW® Technology you'll notice:

- Less fluid in the trocar seals
- Reduced smudging during scope reinsertion
- Fewer interruptions to reclean the scope

For more information, contact your Ethicon Endo-Surgery Sales Professional, call 1-800-USE-ENDO or visit us at ethiconendosurgery.com/xceloptiview



*Compared to most popular competitive 5mm trocars — Applied CQ004, CTR03, CFF03, CTF03, CFH3, CFR03; Covidien NB5STF, VS101005; Taut 40510. Data on file.
©2011 Ethicon Endo-Surgery, Inc. All Rights Reserved. DSL 11-0759.

EES Access

BLADELESS

5mm ENDOPATH® XCEL™ Bladeless Trocars



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
NEW	2B5ST	Stability	75mm	Yes	6
NEW	2B5LT	Stability	100mm	Yes	6
NEW	2B5XT	Stability	150mm	Yes	6
	B5ST	Stability	75mm	No	6
	B5LT	Stability	100mm	No	6
	B5XT	Stability	150mm	No	6
	B5SP*	Smooth	75mm	No	6
	B5LP*	Smooth	100mm	No	6

* Discontinued September 1, 2011

5mm ENDOPATH® XCEL™ Bladeless Trocar with Handle



CODE	SLEEVE	SHAFT LENGTH	QUANTITY
B5LPH*	Smooth	100mm	6

* Discontinued September 1, 2011

8mm ENDOPATH® XCEL™ Bladeless Trocar



CODE	SLEEVE	SHAFT LENGTH	QUANTITY
B8LT	Stability	100mm	6



EES Access

11mm ENDOPATH® XCEL™ Bladeless Trocars



CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
B11LT	Stability	100mm	No	6
B11LP	Smooth	100mm	No	6

11mm ENDOPATH® XCEL™ Bladeless Trocars with Handle



CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
B11LTH	Stability	100mm	No	6
B11LPH	Smooth	100mm	No	6



12mm ENDOPATH® XCEL™ Bladeless Trocars



CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
B12SRT	Stability	75mm	No	6
B12LT	Stability	100mm	No	6
B12XT	Stability	150mm	No	6
B12LP	Smooth	100mm	No	6

EES Access

12mm ENDOPATH® XCEL™ Bladeless Trocars with Handle



CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
B12LTH	Stability	100mm	No	6
B12LPH	Smooth	100mm	No	6

15mm ENDOPATH® XCEL™ Bladeless Trocar



CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
B15LT	Stability	100mm	No	6

EES Access

DILATING TIP

5mm ENDOPATH® XCEL™ Dilating Tip Trocars



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
NEW	2D5ST	Stability	75mm	Yes	6
NEW	2D5LT	Stability	100mm	Yes	6
	D5ST	Stability	75mm	No	6
	D5LT	Stability	100mm	No	6

11mm ENDOPATH® XCEL™ Dilating Tip Trocar



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
	D11LT	Stability	100mm	No	6

12mm ENDOPATH® XCEL™ Dilating Tip Trocars



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
	D12LT	Stability	100mm	No	6
	D12XT	Stability	150mm	No	6

EES Access

SLEEVES

5mm ENDOPATH® XCEL™ Universal Sleeves

Can be used with both XCEL™
Bladeless and Dilating Tip Trocar



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
NEW	2CB5ST	Stability	75mm	Yes	6
NEW	2CB5LT	Stability	100mm	Yes	6
	CB5ST	Stability	75mm	No	6
	CB5LT	Stability	100mm	No	6

11mm ENDOPATH® XCEL™ Universal Sleeve

Can be used with both XCEL™
Bladeless and Dilating Tip Trocar



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
	CB11LT	Stability	100mm	No	6

12mm ENDOPATH® XCEL™ Universal Sleeve

Can be used with both XCEL™
Bladeless and Dilating Tip Trocar



	CODE	SLEEVE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
	CB12LT	Stability	100mm	No	6

EES Access

BLUNT TIP


12mm ENDOPATH® XCEL™ Blunt Tip Trocar



CODE	SHAFT LENGTH	OPTIVIEW® TECHNOLOGY	QUANTITY
H12LP	100mm	No	6

NEEDLES

ENDOPATH® Insufflation Needles Pneumoperitoneum



CODE	SHAFT LENGTH	QUANTITY
PN120	120mm	12
PN150	150mm	12

ENDOPATH® Insufflation Needle Ultra Veress



CODE	SHAFT LENGTH	QUANTITY
UV120	120mm	12

EES Access

MINI/MICRO TROCARS

2mm/3mm ENDOPATH® Mini/Micro Trocars Nonoptical



CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
23NBS	Smooth	2mm/3mm	65mm	6
23NBL	Smooth	2mm/3mm	100mm	6

3mm ENDOPATH® Mini/Micro Trocar



CODE	DESCRIPTION	SHAFT DIAMETER	QUANTITY
AN3MM	3mm Access Needle	3mm	6

5mm ENDOPATH® Nonshielded Trocars



CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
355NS	Nonshielded	5mm	65mm	6
35LNS	Nonshielded	5mm	100mm	6

EES Access



THORACIC TROCARS

10mm/12mm ENDOPATH® Rigid Thoracic Trocar



CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
TT012	Thoracic	10mm/12mm	10mm/12mm	6

FLEXIPATH® Flexible Thoracic Trocars



CODE	OBTURATOR AND SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
FP007	Flexible	7mm	—	6
FP015	Flexible	15mm	—	6
FP020	Flexible	20mm	—	6

FLEXIPATH® Flexible Thoracic Trocar Pack




CODE	DESCRIPTION	QUANTITY
FPK02	Contains: (1) 15mm Obturator and (3) 15mm Flexible Sleeves	6

EES Access

RESPONSABLE TROCARS

ENDOPATH® Bladeless, Obturator and Housing

To be used with ENDOPATH®
Reusable Sleeves

CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
355NA	Resposable	5mm	75mm	6
35LNA	Resposable	5mm	100mm	6
511NA	Resposable	10mm/11mm	100mm	6
512NA	Resposable	10mm/12mm	100mm	6

ENDOPATH® Dilating Tip, Obturator and Housing

To be used with ENDOPATH®
Reusable Sleeves

CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
355DA	Resposable	5mm	75mm	6
35LDA	Resposable	5mm	100mm	6
511DA	Resposable	10mm/11mm	100mm	6
512DA	Resposable	10mm/12mm	100mm	6

ENDOPATH® Blunt Tip, Obturator, Housing and Olive Plug



CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
512BA	Resposable	10mm/12mm	100mm	6

EES Access

HOUSING UNITS AND SLEEVES

ENDOPATH® Disposable Housing Units



CODE	SHAFT DIAMETER	QUANTITY
355HR	5mm	6
511HR	10mm/11mm	6
512HR	10mm/12mm	6

ENDOPATH® Reusable Sleeves



CODE	SLEEVE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
355RT	Stability	5mm	75mm	1
351RT	Stability	5mm	100mm	1
511RT	Stability	10mm/11mm	100mm	1
512RT	Stability	10mm/12mm	100mm	1

SINGLE-SITE LAPAROSCOPY

Single-Site Laparoscopy



CODE	DESCRIPTION	QUANTITY
SSL01FLR	Single-Site Laparoscopy Fixed-Length Retractor (Small)	6
SSL02FLR	Single-Site Laparoscopy Fixed-Length Retractor (Medium)	6
SSL3P	Single-Site Laparoscopy Seal Cap with Accessories	3


EES Access

HAND-ASSISTED LAPAROSCOPY (HALS)

ENDOPATH® DEXTRUS™ HALS Retractors



FLR01



FLR02



FLR03


CODE	DESCRIPTION	ABDOMINAL WALL THICKNESS	QUANTITY
FLR01	Fixed-Length Access Retractor (Small)	<4cm	6
FLR02	Fixed-Length Access Retractor (Medium)	4-7cm	6
FLR03	Fixed-Length Access Retractor (Large)	>7cm	6

ENDOPATH® DEXTRUS™ HALS Products



CODE	DESCRIPTION	QUANTITY
HAP02	Seal Cap Assembly with Accessories	3






Stapling

How much unintentional tip movement does your endocutter have when firing?

Introducing the gold standard in stability: The NEW ECHELON FLEX™ Powered ENDOPATH® Stapler

Allows the surgeon to reduce movement at the end effector by at least 63% during transection of thick tissue*

- **Reduces unwanted movement** of the distal tip to protect targeted tissue and surrounding fragile organs and tissue from damage
- **Results in smooth, consistent firing** even in thick tissue and at awkward angles
- **Reduced force to fire** leads to less movement of the tip on vital structures, potentially resulting in less trauma



See for yourself what a difference the new Powered Echelon Flex™ Stapler can make in your procedures.
Visit us at ethiconendosurgery.com/poweredechelexonflex

*Benchtop testing on porcine stomach (2.5mm to 4mm thick). Surgeons (n=19) fired each instrument/reload once: PSE60A/ECR60G, 030449/t30459 and EGIAUSTND/EGIA60AMT. Distal tip motion measurement during the firing cycle showed a 63% reduction in tip movement of PSE60A/ECR60G vs the other two devices.
© 2011 Ethicon Endo-Surgery, Inc. All Rights Reserved. DSL 11-0865.

POWERED ENDOCUTTERS

Powered ECHELON FLEX™ 45mm Stapler

COMING SOON






CODE	ARTICULATION	SHAFT LENGTH	QUANTITY
PCE45A	Yes	280mm	3
PSE45A	Yes	340mm	3
PLE45A	Yes	440mm	3

EES Stapling


ECHELON™ 45mm Reloads

For use with all ECHELON™ 45mm Endoscopic Staplers
(SC45A, EC45A, EC45AL, SC45, EC45, ECLG45)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
ECR45M	Gray	6	Mesentery	2.0mm	.75mm	12
ECR45W	White	6	Vascular	2.5mm	1.0mm	12
ECR45B	Blue	6	Standard	3.5mm	1.5mm	12
ECR45D	Gold	6	Standard/Thick	3.8mm	1.8mm	12
ECR45G	Green	6	Thick	4.1mm	2.0mm	12



Powered ECHELON FLEX™ 60mm Stapler



CODE	ARTICULATION	SHAFT LENGTH	QUANTITY
PCE60A	Yes	280mm	3
PSE60A	Yes	340mm	3
PLE60A	Yes	440mm	3

ECHELON™ 60mm Reloads

For use with all ECHELON™ 60mm Endoscopic Staplers
(SC60A, EC60A, LONG60A, SC60, EC60, LONG60)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
ECR60M	Gray	6	Mesentery	2.0mm	.75mm	12
ECR60W	White	6	Vascular	2.5mm	1.0mm	12
ECR60B	Blue	6	Standard	3.5mm	1.5mm	12
ECR60D	Gold	6	Standard/Thick	3.8mm	1.8mm	12
ECR60G	Green	6	Thick	4.1mm	2.0mm	12

EES Stapling

ECHELON FLEX™ 45mm Endocutters Reloads Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
SC45A	Yes	280mm	12	3
EC45A	Yes	340mm	12	3
EC45AL	Yes	440mm	12	3


ECHELON™ 45mm Endocutters Reloads Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
SC45	No	280mm	12	3
EC45	No	340mm	12	3
ECLG45	No	440mm	12	3

ECHELON™ 45mm Reloads


For use with all ECHELON™ 45mm Endoscopic Staplers
(SC45A, EC45A, EC45AL, SC45, EC45, ECLG45)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
ECR45M	● Gray	6	Mesentery	2.0mm	.75mm	12
ECR45W	○ White	6	Vascular	2.5mm	1.0mm	12
ECR45B	● Blue	6	Standard	3.5mm	1.5mm	12
ECR45D	○ Gold	6	Standard/Thick	3.8mm	1.8mm	12
ECR45G	● Green	6	Thick	4.1mm	2.0mm	12


EES Stapling

ECHELON FLEX™ 60mm Endocutters Cartridges Sold Separately








CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
SC60A	Yes	280mm	12	3
EC60A	Yes	340mm	12	3
LONG60A	Yes	440mm	12	3

ECHELON™ 60mm Endocutters Cartridges Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
SC60	No	280mm	12	3
EC60	No	340mm	12	3
LONG60	No	440mm	12	3

ECHELON™ 60mm Reloads For use with all ECHELON™ 60mm Endoscopic Staplers (SC60A, EC60A, LONG60A, SC60, EC60, LONG60)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
ECR60M	● Gray	6	Mesentery	2.0mm	.75mm	12
ECR60W	○ White	6	Vascular	2.5mm	1.0mm	12
ECR60B	● Blue	6	Standard	3.5mm	1.5mm	12
ECR60D	○ Gold	6	Standard/Thick	3.8mm	1.8mm	12
ECR60G	● Green	6	Thick	4.1mm	2.0mm	12

EES Stapling

ETS Flex 45mm Endocutters Reloads Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
CTS45	Yes	240mm	8	3
ATS45	Yes	340mm	8	3
LONG45A	Yes	440mm	8	3

ETS Flex 45mm No-Knife Endocutters Reloads Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
CTS45NK	Yes	240mm	8	3
ATS45NK	Yes	340mm	8	3

ETS 45mm Endocutter Reloads Sold Separately



CODE	ARTICULATION	SHAFT LENGTH	MAXIMUM FIRINGS	QUANTITY
ETS45	No	340mm	8	3


ETS 45mm Reloads For use with all ETS 45mm Endoscopic Linear Cutters (CTS45, ATS45, LONG45A, CTS45NK, ATS45NK, ETS45)





CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
6R45M	Gray	6	Mesentery	2.0mm	.85mm	12
TR45W	White	6	Vascular	2.5mm	1.0mm	12
6R45B	Blue	6	Standard	3.5mm	1.5mm	12
TR45G	Green	4	Standard/Thick	4.1mm	2.0mm	12



EES Stapling



EZ 45mm Endocutters



CODE	COLOR	ARTICULATION	SHAFT LENGTH	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
EZ45B	 Blue	No	300mm	ZR45B	8	3
EZ45G	 Green	No	300mm	ZR45G	8	3


EZ 45mm Reloads For use with EZ 45mm Endoscopic Linear Cutters (EZ45B, EZ45G)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	MAXIMUM FIRINGS	QUANTITY
ZR45B	 Blue	4	Standard	3.5mm	1.5mm	8	12
ZR45G	 Green	4	Thick	4.5mm	2.0mm	8	12


EES Stapling

ETS Flex 35mm Endocutters ETS Flex 35mm Articulating Endoscopic Linear Cutters



CODE	ARTICULATION	SHAFT LENGTH	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
ATW35	Yes	340mm	TR35W	8	3
ATB35	Yes	340mm	TR35B	8	3

ETS 35mm Straight Endocutters



CODE	ARTICULATION	SHAFT LENGTH	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TSW35	No	340mm	TR35W	8	3
TSB35	No	340mm	TR35B	8	3

ETS 35mm Reloads For use with ETS 35mm Endoscopic Linear Cutters (ATW35, ATB35, TSW35, TSB35)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
TR35W	○ White	6	Vascular	2.5mm	1.0mm	12
TR35B	● Blue	6	Standard	3.5mm	1.5mm	12

EES Stapling


LINEAR CUTTERS

55mm Linear Cutter, Selectable Staple Height Reloads Sold Separately



CODE	SELECTABLE COLOR	TISSUE TYPE	MAXIMUM FIRINGS	QUANTITY
NTLC55	● ○ ●	Selectable	12	3

55mm Linear Cutter, Selectable Reload



CODE	SELECTABLE COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
SR55	●	6	Standard	4.3mm	1.5mm	12
	○	6	Standard/Thick	4.3mm	1.8mm	12
	●	6	Thick	4.3mm	2.0mm	12

75mm Linear Cutter, Selectable Staple Height Reloads Sold Separately



CODE	SELECTABLE COLOR	TISSUE TYPE	MAXIMUM FIRINGS	QUANTITY
NTLC75	● ○ ●	Selectable	12	3

75mm Linear Cutter, Selectable Reload



CODE	SELECTABLE COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
SR75	●	6	Standard	4.3mm	1.5mm	12
	○	6	Standard/Thick	4.3mm	1.8mm	12
	●	6	Thick	4.3mm	2.0mm	12

EES Stapling

55mm Linear Cutters PROXIMATE® 55mm Linear Cutters with Safety Lockout



TVC55





TLC55



TCT55

CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TVC55		White Vascular	2.5mm	1.0mm	TVR55	8	3
TLC55		Blue Standard	3.85mm	1.5mm	TCR55	8	3
TCT55		Green Thick	4.5mm	2.0mm	TRT55	8	3

55mm Linear Cutter Reloads For use with all 55mm Linear Cutters (TVC55, TLC55, TCT55)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
TVR55		4	Vascular	2.5mm	1.0mm	12
TCR55		4	Standard	3.85mm	1.5mm	3
TRT55		4	Thick	4.50mm	2.0mm	3

75mm Linear Cutters PROXIMATE® 75mm Linear Cutters with Safety Lockout



TLC75






TCD75



TCT75

CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLC75		Blue Standard	3.85mm	1.5mm	TCR75	8	3
TCD75		Gold Standard/Thick	4.2mm	1.8mm	TRD75	8	3
TCT75		Green Thick	4.5mm	2.0mm	TRT75	8	3

75mm Linear Cutter Reloads For use with all 75mm Linear Cutters (TLC75, TCD75, TCT75)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
TCR75		4	Standard	3.85mm	1.5mm	12
TRD75		4	Standard/Thick	4.2mm	1.8mm	12
TRT75		4	Thick	4.5mm	2.0mm	12

EES Stapling

100mm Linear Cutters PROXIMATE® 100mm Linear Cutters with Safety Lockout





TLC10



TCT10

CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLC10	● Blue	Standard	3.85mm	1.5mm	TCR10	8	3
TCT10	● Green	Thick	4.5mm	2.0mm	TRT10	8	3

100mm Linear Cutter Reloads For use with all 100mm Linear Cutters (TLC10, TCT10)



CODE	COLOR	ROWS	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
TCR10	● Blue	4	Standard	3.85mm	1.5mm	12
TRT10	● Green	4	Thick	4.5mm	2.0mm	12

EES Stapling

CIRCULAR STAPLERS

Endoscopic Curved Intraluminal Staplers



ECS21A



ECS25A



ECS29A



ECS33A

CODE	LUMEN	SHAFT LENGTH	HEAD DIAMETER	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
ECS21A	12.4mm	28cm	21mm	5.5mm	1.0-2.5mm Adjustable	3
ECS25A	16.4mm	28cm	25mm	5.5mm	1.0-2.5mm Adjustable	3
ECS29A	20.4mm	28cm	29mm	5.5mm	1.0-2.5mm Adjustable	3
ECS33A	24.4mm	28cm	33mm	5.5mm	1.0-2.5mm Adjustable	3

Curved Intraluminal Staplers



CDH21A



CDH25A




CDH29A



CDH33A

CODE	LUMEN	SHAFT LENGTH	HEAD DIAMETER	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
CDH21A	12.4mm	18cm	21mm	5.5mm	1.0-2.5mm Adjustable	3
CDH25A	16.4mm	18cm	25mm	5.5mm	1.0-2.5mm Adjustable	3
CDH29A	20.4mm	18cm	29mm	5.5mm	1.0-2.5mm Adjustable	3
CDH33A	24.4mm	18cm	33mm	5.5mm	1.0-2.5mm Adjustable	3

Straight Intraluminal Staplers



SDH21A



SDH25A



SDH29A




SDH33A

CODE	LUMEN	SHAFT LENGTH	HEAD DIAMETER	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
SDH21A	12.4mm	18cm	21mm	5.5mm	1.0-2.5mm Adjustable	3
SDH25A	16.4mm	18cm	25mm	5.5mm	1.0-2.5mm Adjustable	3
SDH29A	20.4mm	18cm	29mm	5.5mm	1.0-2.5mm Adjustable	3
SDH33A	24.4mm	18cm	33mm	5.5mm	1.0-2.5mm Adjustable	3

EES Stapling

ENDOPATH® ILS Circular Sizer Set



CODE	DESCRIPTION	QUANTITY
CSS	Bowel Sizers for the ENDOPATH® ILS Curved Staplers, the set includes one each of a 25mm, 29mm and 33mm sizer	1

PROCEDURE STAPLERS

PPH Procedure Set




CODE	DESCRIPTION	QUANTITY
PPH03	Includes Hemorrhoidal Circular Stapler 33mm, suture threader, circular anal dilator and purse string anoscope	3

STARR Procedure Set



CODE	DESCRIPTION	QUANTITY
STR10	For Stapled Transanal Rectal Resection (STARR)	3

Reusable Purse String Clamp




CODE	DESCRIPTION	QUANTITY
EH40	—	1

EES Stapling

LINEAR STAPLERS

TX 30mm Linear Stapler PROXIMATE™ TX 30mm Reloadable Linear Stapler



CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	QUANTITY
TX30V	○ White	Vascular	2.5mm	1.0mm	XR30V	3
XR30V	○ White	Vascular	2.5mm	1.0mm	—	12

TX 30mm Linear Staplers PROXIMATE™ TX 30mm Reloadable Linear Staplers



CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	QUANTITY
TX30B	● Blue	Standard	3.5mm	1.5mm	XR30B	3
TX30G	● Green	Thick	4.8mm	2.0mm	XR30G	3


TX 30mm Linear Stapler Reloads For use with PROXIMATE™ TX 30mm Reloadable Linear Staplers (TX30B, TX30G)





CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
XR30B	● Blue	Standard	3.5mm	1.5mm	12
XR30G	● Green	Thick	4.8mm	2.0mm	12

EES Stapling



TX 60mm Linear Staplers PROXIMATE® TX 60mm Reloadable Linear Staplers



CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	QUANTITY
TX60B	 Blue	Standard	3.5mm	1.5mm	XR60B	3
TX60G	 Green	Thick	4.8mm	2.0mm	XR60G	3

TX 60mm Linear Stapler Reloads For use with PROXIMATE® TX 60mm Reloadable Linear Staplers (TX60B, X60G)



CODE	COLOR	TISSUE TYPE	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
XR60B	 Blue	Standard	3.5mm	1.5mm	12
XR60G	 Green	Thick	4.8mm	2.0mm	12

EES Stapling

TL 30mm Linear Stapler with Dial

TL Reloadable Linear Stapler with Adjustable Staple Height and Reload



CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TL30	30mm TL Reloadable Linear Stapler	1.0-2.5mm	TR30	4	3
TR30	—	1.0-2.5mm	—	—	12

TL 60mm Linear Stapler with Dial


TL Reloadable Linear Stapler with Adjustable Staple Height and Reload



CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TL60	60mm TL Reloadable Linear Stapler	1.0-2.5mm	TR60	4	3
TR60	—	1.0-2.5mm	—	—	12

TL 90mm Linear Stapler with Dial

TL Reloadable Linear Stapler, Heavy Wire with Adjustable Staple Height and Reload



CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TL90	90mm TL Reloadable Linear Stapler	1.0-2.5mm	TR90	4	3
TR90	—	1.0-2.5mm	—	—	12

EES Stapling

TLH 30mm Linear Stapler with Dial TLH Reloadable Linear Stapler with Adjustable Staple Height




CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLV30	30mm TLH Reloadable Linear Stapler	1.0mm	TRV30	4	3
TRV30	For use with TLV30	1.0mm	—	—	12

TLH 30mm Linear Stapler with Dial TLH Reloadable Linear Stapler, Heavy Wire with Adjustable Staple Height




CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLH30	30mm TLH Reloadable Linear Stapler	1.0-2.5mm	TRH30	4	3
TRH30	For use with TLH30	1.0-2.5mm	—	—	12

TLH 60mm Linear Stapler with Dial TLH Reloadable Linear Stapler, Heavy Wire with Adjustable Staple Height



CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLH60	60mm TLH Reloadable Linear Stapler	1.5-2.5mm	TRH60	4	3
TRH60	For use with TLH60	1.5-2.5mm	—	—	12

TLH 90mm Linear Stapler with Dial TLH Reloadable Linear Stapler, Heavy Wire with Adjustable Staple Height



CODE	DESCRIPTION	STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
TLH90	90mm TLH Reloadable Linear Stapler	1.5-2.5mm	TRH90	4	3
TRH90	For use with TLH90	1.5-2.5mm	—	—	12

EES Stapling


SKIN STAPLERS

Rotating-Head Skin Staplers PROXIMATE® Rotating-Head Skin Staplers with Stainless Steel Staples




CODE	SIZE	STAPLES	QUANTITY
PRR35	Regular	35	6
PRW35	Wide	35	6

Fixed-Head Skin Staplers PROXIMATE® Plus MD Fixed-Head Skin Staplers with Multidirectional Release



CODE	SIZE	STAPLES	QUANTITY
PMR35	Regular	35	6
PMW35	Wide	35	6

Fixed-Head Skin Staplers PROXIMATE® PX Fixed-Head Skin Staplers with Multidirectional Release, Pistol Grip



CODE	SIZE	STAPLES	QUANTITY
PXR35	Regular	35	6
PXW35	Wide	35	6

Staple Extractor






CODE	SIZE	QUANTITY
PSX	Squeeze Handle Staple Extractor	12

EES Stapling

CURVED CUTTER STAPLERS

Contour® 40mm Linear Staplers



CODE	COLOR	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	MAXIMUM FIRINGS	QUANTITY
CS40B	Blue	3.5mm	1.44mm	CR40B	6	3
CS40G	Green	4.7mm	2.0mm	CR40G	6	3

Contour® 40mm Linear Stapler Reloads For use with (CS40B, CS40G)



CODE	COLOR	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	QUANTITY
CR40B	Blue	3.5mm	1.44mm	6
CR40G	Green	4.7mm	2.0mm	6

55mm Articulating Linear Staplers Single-Patient Firings



CODE	COLOR	OPEN STAPLE HEIGHT	CLOSED STAPLE HEIGHT	RELOAD INCLUDED	QUANTITY
AX55B	Blue	3.5mm	1.5mm	Not Reloadable	3
AX55G	Green	4.8mm	2.0mm	Not Reloadable	3

EES Stapling

Reach Deeper, See Better¹

Contour® Curved Cutter Staplers.

- Features **simultaneous stapling and cutting**
- Delivers **four curvilinear rows of staples with a single cut**.
- Unique curved-head design lets you **access the lower pelvis and enhances your visibility** of the targeted area.
- A single device may be used to **fire up to six times** in a single procedure
- Offers an indication for **skeletal muscle**, including steps to use for **muscle stapling** in below- and above-knee amputations and sarcoma cases.



¹Unique curved-head design gives lower pelvic access than 30mm linear staplers (PISC) and enhanced visibility (of anatomic structures). Data on file at EES.

For more information, contact your
Ethicon Endo-Surgery Sales Professional or call 1-800-USE-ENDO.







Ligation

EES Ligation



ENDOSCOPIC MULTIPLE-CLIP APPLIERS

LIGAMAX™ 5 Rotating Multiple-Clip Applier





CODE	DIAMETER	CLIP SIZE	NUMBER OF CLIPS	QUANTITY
EL5ML	5mm	Medium/Large	15	3

LIGACLIP® Rotating Multiple-Clip Applier



CODE	DIAMETER	CLIP SIZE	NUMBER OF CLIPS	QUANTITY
ER320	10mm	Medium/Large	20	3


LIGACLIP® Rotating Multiple-Clip Applier



CODE	DIAMETER	CLIP SIZE	NUMBER OF CLIPS	QUANTITY
ER420	12mm	Large	20	3

ENDOSCOPIC SINGLE-CLIP APPLIERS

LIGACLIP® Multi-Patient-Use Single-Clip Appliers (Uses LIGACLIP® EXTRA Stainless Steel or Titanium Ligating Clips)



CODE	CLIP SIZE	QUANTITY
EL214	Medium	1
EL314	Medium/Large	1
EL414	Large	1

EES Ligation

OPEN MULTIPLE-CLIP APPLIERS

LIGACLIP® Multiple-Clip Appliers Single Use, Not Reloadable



MCS20



MSM20



MCM20



MCM30




MCL20

CODE	LENGTH	CLIP SIZE	NUMBER OF CLIPS	QUANTITY
MCS20	23.8cm	Small	20	6
MSM20	23.8cm	Medium	20	6
MCM20	29.2cm	Medium	20	6
MCM30	29.2cm	Medium	30	6
MCL20	33.7cm	Large	20	6

The enhanced Ligamax™ 5. Superior clip security*

5mm wide. 10mm good.

- **Delivers the same security** of medium/large clips in a small clip applicator
 - Your patients will appreciate all of the benefits of a **smaller incision**
- That's why it's the market-leading 5mm clip applicator



So what's stopping you from switching?
Visit us at ethiconendosurgery.com/tryligamax5

*Only on file. Benchtop testing using elastomers. Testing comparing JIGAMAX™ and ENDO CLIP™ II (Ethicon Endo-Surgery, Inc.) and ENDO CLIP™ II (Ethicon Endo-Surgery, Inc.). All Rights Reserved. DS1-11-0002

EES Ligation

OPEN SINGLE-CLIP APPLIERS


LIGACLIP® Multi-Patient-Use Single-Clip Appliers — Blue




LX105



LX107



LX110



LX130

CODE	DESCRIPTION	HANDLE COLOR	LENGTH	CLIP SIZE	QUANTITY
LX105	—	Blue	14.6cm	Small	1
LX107	—	Blue	19.1cm	Small	1
LX110	—	Blue	26.7cm	Small	1
LX130	Right Angle	Blue	26.7cm	Small	1

LIGACLIP® Multi-Patient-Use Single-Clip Appliers — Silver



LX205



LX207



LX210



LX220



LX230

CODE	DESCRIPTION	HANDLE COLOR	LENGTH	CLIP SIZE	QUANTITY
LX205	—	Silver	14.6cm	Medium	1
LX207	—	Silver	19.1cm	Medium	1
LX210	—	Silver	26.7cm	Medium	1
LX220	Angled Jaw	Silver	26.7cm	Medium	1
LX230	Right Angle	Silver	26.7cm	Medium	1

LIGACLIP® Multi-Patient-Use Single-Clip Appliers — Green



LC307



LC310



LC320



LC330


CODE	DESCRIPTION	HANDLE COLOR	LENGTH	CLIP SIZE	QUANTITY
LC307	—	Green	19.1cm	Medium/Large	1
LC310	—	Green	26.7cm	Medium/Large	1
LC320	Angled Jaw	Green	26.7cm	Medium/Large	1
LC330	Right Angle	Green	26.7cm	Medium/Large	1

EES Ligation


LIGACLIP® Multi-Patient-Use Single-Clip Appliers — Yellow



LC407



LC410



LC420



LC430

CODE	DESCRIPTION	HANDLE COLOR	LENGTH	CLIP SIZE	QUANTITY
LC407	—	Yellow	19.1cm	Large	1
LC410	—	Yellow	26.7cm	Large	1
LC420	Angled Jaw	Yellow	26.7cm	Large	1
LC430	Right Angle	Yellow	26.7cm	Large	1

LIGACLIP® Ligating Clip Cartridge Base



CODE	DESCRIPTION	QUANTITY
LC800	LIGACLIP® Ligating Clip Cartridge Base	1

SINGLE-CLIPS

LIGACLIP® EXTRA Titanium Ligating Clips

For use with both Endoscopic and Open/Conventional Single-Clip Reusable Appliers



LT100



LT102



LT200



LT202



LT300



LT400

CODE	COLOR	CLIP SIZE	NUMBER OF CLIPS	QUANTITY
LT100	Blue	Small	6	36
LT102	Blue	Small	20	15
LT200	Silver	Medium	6	36
LT202	Silver	Medium	20	15
LT300	Green	Medium/Large	6	18
LT400	Yellow	Large	6	18





Endo Devices

EES Endo Devices

POUCH


ENDOPOUCH® Specimen Retrieval Bags




CODE	DESCRIPTION	DIAMETER	QUANTITY
POUCH	4"x 6" ENDOPOUCH® Specimen Retrieval Bag	10mm	6

ELECTROSURGICAL IRRIGATION

ENDOPATH® Laparoscopic Bipolar Forceps



EBF01



EBF02

CODE	JAW	DIAMETER	QUANTITY
EBF01	Macro	5mm	5
EBF02	Micro	5mm	5

ENDOPATH® PROBE PLUS® II Handles



EPH01



EPH02



EPH03



EPH04

CODE	CONTROL	GRIP	QUANTITY
EPH01	Foot	Pistol	6
EPH02	Hand	Pistol	6
EPH03	Foot	Pencil	6
EPH04	Hand	Pencil	6

EES Endo Devices

ENDOPATH® PROBE PLUS® II Hollow Tip Electrodes



CODE	ELECTRODE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
EPS01	Hook	5mm	34cm	6
EPS02	Spatula	5mm	34cm	6
EPS03	Right Angle	5mm	34cm	6
EPS04	Curved Dissector	5mm	34cm	6

ENDOPATH® PROBE PLUS® II Hollow Tip Electrodes



CODE	ELECTRODE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
EPS05	Hook	5mm	29cm	6
EPS06	Spatula	5mm	29cm	6
EPS07	Right Angle	5mm	29cm	6
EPS08	Curved Dissector	5mm	29cm	6

ENDOPATH® PROBE PLUS® II Irrigation



CODE	SUCTION	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
EPS10	Pool/Sump	10mm	34cm	6
EPS11	Open End	10mm	34cm	6

EES Endo Devices

5mm ENDOPATH® PROBE PLUS® II Canulas



CODE	TYPE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
EP513	Flexible Fiber	5mm	29cm	6

10mm ENDOPATH® PROBE PLUS® II Canulas



CODE	TYPE	SHAFT DIAMETER	SHAFT LENGTH	QUANTITY
EP512	Stone Retrieval	10mm	34cm	6

ENDOPATH® PROBE PLUS® II Tubing



CODE	DESCRIPTION	QUANTITY
EPT01	PROBE PLUS® II, Dual-Bottle Irrigation Tube Set	6
EPT03	PROBE PLUS® II, Pressurized Bag Irrigation Tube Set with Spike	6

ENDOPATH® Active Cords



CODE	DESCRIPTION	PRONGS	QUANTITY
EBC02	Reusable Coaxial Pin Connector Bipolar Cord	1	1
EBC05	Active Return Cord for Bipolar Instruments, Fixed Pin Connector	2	1

EES Endo Devices


ENDO INSTRUMENTS

ENDOPATH® 5mm Instruments




CODE	DESCRIPTION	QUANTITY
5BB	Babcock, 5mm	6
5DCD	Dissector, Curved, with Monopolar Cautery, 5mm	6
5DCS	Scissors, Curved, with Monopolar Cautery, 5mm	6
5DSG	Grasper, 5mm	6

ENDOPATH® 5mm Instruments — Blunt Tip



CODE	DESCRIPTION	QUANTITY
BTDO5	Dissector, Blunt Tip, 5mm, 3 per Pouch, 12 Pouches per Unit	36


ENDOPATH® 10mm Instruments



CODE	DESCRIPTION	QUANTITY
10AG	Anvil Grasper, 10mm	6
10BB	Babcock, 10mm	6

EES Endo Devices


ENDOPATH® 10mm Instruments — Blunt Tip



CODE	DESCRIPTION	QUANTITY
BCD10	Dissector, Blunt Cherry, 10mm, 3 per Pouch, 12 Pouches per Unit	36


ENDOSUTURE

LAPRA-TY® Absorbable Suture Products




CODE	DESCRIPTION	SHAFT	QUANTITY
KA200	LAPRA-TY® Absorbable Suture Clip Applier, Reusable, for use with 11mm or 12mm Diameter Trocar	33cm	1
NEW KA200L	Long LAPRA-TY® Absorbable Suture Clip Applier, Reusable, for use with 11mm or 12mm Diameter Trocar	45cm	1

LAPRA-TY® Absorbable Suture Products



CODE	DESCRIPTION	QUANTITY
XC200	LAPRA-TY® Absorbable Suture Clips for use with 2-0, 3-0 and 4-0 Coated VICRYL® Suture, 6 Clips per Cartridge	6


Suture Assistant



CODE	DESCRIPTION	QUANTITY
SW100	Endoscopic Suturing Device	4

EES Endo Devices

Suture Assistant



SW110



SW112



SW120



SW122



SW222

CODE	DESCRIPTION	QUANTITY
SW110	Suture Reload Unit with 6" ETHIBOND™ EXCEL (size 0) and EEN (CANOE) Needle	10
SW112	Suture Reload Unit with 6" ETHIBOND™ EXCEL (size 2-0) and EEN (CANOE) Needle	10
SW120	Suture Reload Unit with 6" ETHIBOND™ EXCEL (size 0) and SH Needle	10
SW122	Suture Reload Unit with 6" ETHIBOND™ EXCEL (size 2-0) and SH Needle	10
SW222	Suture Reload Unit with 6" VICRYL® (size 2-0) and SH Needle	10

ENDOLOOP® Ligature




EJ10G



EZ10G

CODE	DESCRIPTION	QUANTITY
EJ10G	ENDOLOOP® Coated VICRYL® Ligature, 0, 18"	12
EZ10G	ENDOLOOP® PDS II Ligature, 0, 18"	12


ENDOKNOT® Suture



EX10G



JK10G



ZK10G


CODE	DESCRIPTION	QUANTITY
EX10G	ETHIBOND™ EXCEL Suture, ST-3 Needle, 0, 42"	12
JK10G	Coated, VICRYL® Suture, ST-3 Needle, 0, 42"	12
ZK10G	PDS II Suture, ST-3 Needle, 0, 42"	12

EES Endo Devices

Needle/Suture Combinations



ESS15



EC11




EC12

CODE	DESCRIPTION	QUANTITY
ESS15	ENDOHOLDER for use with 5mm Diameter ENDOPATH [®] Trocar, Reusable	1
EC11	ETHIBOND [™] EXCEL Suture, EN-3 needle, 0, 44"	12
EC12	ETHIBOND [™] EXCEL Suture, EN-3 needle, 2-0, 44"	12


Port Closure Needles



J00JG



JTN0G



Z00JG

CODE	DESCRIPTION	QUANTITY
J00JG	Port Closure Needle, Coated VICRYL [®] J needle, 0, 18"	12
JTN0G	Port Closure Needle, Coated VICRYL [®] TN needle, 0, 27"	12
Z00JG	Port Closure Needle, PDS II, J, 0, 18"	12

MISCELLANEOUS

Hernia Fixation



CODE	DESCRIPTION	QUANTITY
EMS	10mm Hernia Tacker	6

EES Endo Devices

Laparoscopic Needle Holders




CODE	DESCRIPTION	QUANTITY
E705R	ENDOPATH® Needle Holder, Reusable, 5mm, 30cm Length	1
E705RL	ENDOPATH® Needle Holder, Reusable, 5mm, 45cm Length	1



CODE	DESCRIPTION	QUANTITY
SRNH1	Self-righting Needle Holder, Reusable, 5mm, 30cm Length	1
SRNH1L	Self-righting Needle Holder, Reusable, 5mm, 45cm Length	1

Uterine Manipulators



CODE	DESCRIPTION	QUANTITY
UM201*	ENDOPATH® Uterine Manipulator	5
UM202*	ENDOPATH® Uterine Manipulator, Long	5

* Discontinued December 31, 2011



CATEGORIES

- Sutures: Absorbable
- Sutures: Non-Absorbable
- Sutures: By Needle Type
- Adhesives
- Drains
- Hernia Repair
- ETHICON SECURESTRAP™
- ETHICON PHYSIOMESH™
- PROCEED™ Mesh
- PVP™ Device
- ULTRAPRO Mesh
- ULTRAPRO Plug
- ULTRAPRO Hernia System
- PROLENE Hernia System
- PROLENE Mesh & PROLENE Soft Mesh
- PROLENE 3D Patch Mesh
- VICRYL Woven Mesh
- MERSILENE Mesh
- FlexHD® Acellular Hydrated Dermis
- Miscellaneous

SPECIFY SEARCH CRITERIA

GO

Advanced Suture Search


Print this Page | Product Catalog in PDF

HERNIA REPAIR: ETHICON SECURESTRAP™

ETHICON SECURESTRAP™ Fixation Device

The ETHICON SECURESTRAP™ Fixation Device is intended for fixation of prosthetic material to soft tissues in various minimally invasive and open surgical procedures, such as hernia repairs.

For additional information about Hernia Repair please go to herniasolutions.com



Product Groups

- ETHICON SECURESTRAP™ 5mm Absorbable Strap Fixation Device

ETHICON SECURESTRAP™ 5mm Absorbable Strap Fixation Device

STRAP12	5mm absorbable Fixation Device	6 devices per box
STRAP25	5mm single-use mesh fixation device w/ 25 absorbable straps	6 devices per box

Back to top

[Legal Notice](#) [Privacy Policy](#) [Site Map](#)

By using this site, you agree to our [Legal Notice](#) and [Privacy Policy](#) © Ethicon, Inc. 2009 – Last updated on October 01, 2009

*Trademark

This site is published by Ethicon, Inc. which is solely responsible for its contents. It is intended for healthcare professionals in the United States.

Home » Products » Ethicon Products » ETHICON SECURESTRAP™ Fixation Device

Prescribing Information

About the
Product

Clinical
Information

Related
Information

ETHICON

SECURESTRAP™

5 mm Absorbable Strap Fixation Device

FOR LAPAROSCOPIC HERNIA REPAIR

Low Strap Indicator




Orientation and Positioning System




SECURESTRAP™

Spring-loaded deployment system
with lock-out handle

25 Synthetic Absorbable Dyed Straps



Unique, absorbable "strap" design uses 2 points of fixation to straddle mesh pores and fibers, including large-pore constructions.⁴



Straddling a mesh pore.



Straddling a mesh fiber.

- Delivers secure mesh fixation during the tissue integration period.
- Helps prevent mesh migration.
- Strap design securely fixates a variety of meshes, including large-pore constructions.
- Unique polymer blend is essentially absorbed within 12-18 months.

References:

- 1-2. Data on file, Ethicon, Inc.
3. Ethicon Securestrap 5mm Absorbable Strap Fixation Device: Instructions for Use.
4. Data on file, Ethicon, Inc.

For more information about the ETHICON SECURESTRAP™ Fixation Device and other ETHICON Hernia Solutions products, contact your ETHICON PRODUCTS account representative at 1-800-4ETHICON.

References:

1. Data on file, Ethicon, Inc.

Home » Products » ETHICON Products » ETHICON SECURESTRAP™ Fixation Device


ETHICON SECURESTRAP™ Fixation Device

Prescribing Information


About the Product Clinical Information Professional Resources Related Information

UNIQUELY DESIGNED TO DELIVER CONSISTENT AND SECURE MESH FIXATION


Delivers secure mesh fixation during the tissue integration period³



Unique strap design allows for appropriate tissue capture and penetration at various deployment angles to provide excellent holding strength.²




Uses 2 points of fixation to straddle mesh pores and fibers, securely fixing a variety of meshes including large-pore constructions.²



Low profile strap minimizes foreign material exposed to viscera.³ Absorbable leaving less permanently implanted material in the patient.⁶

ETHICON SECURESTRAP™ video



References:

1. 510K Testing for ORION Mesh Fixation Strength (AST-2009-0312).
2. Angle Fire Competitive Testing (AST- 2010-0199)
3. 28-Day Mesh Fixation Study of the ORION Device to Evaluate Mesh Migration and Tissue Response Using a Swine Model (PSE Accession Number 09-0132).
4. ETHICON SECURESTRAP™ Design Reliability Testing.
5. ETHICON SECURESTRAP™ 5mm Absorbable Strap Fixation Device: Instructions for Use.
6. Prediction of In Vivo Total Absorption Time for ORION Implant Using In Vitro/In Vivo Correlation Degradation Model (AST-2009-0332).
7. Divilio LT: Surgical adhesion development and prevention. Int Surg (2005) 22:2021-2024.

